

# Bioplus

Changing the World of Pipetting



OEM AVAILABLE

## Every Drop of Precision

- New ergonomic design
- Instant volume adjustments
- Separate tip ejector
- Clear digital window on the hand grip



**BIOSYSTEM**<sup>®</sup>  
'ACCURACY & PRECISION'

CE ISO 9001 : 2000 CERTIFICATE CO.

# Bioplus

Biosystem is a compact device with separate tip ejector. Bioplus pipettes enjoys excellent reputation with medical researchers and professional people. It comes with superior feature to achieve highest efficiency. The unit is very user and strain free.



## Features of Bioplus

- New ergonomic design
- It reduced probability of repetitive strain injure
- Instant volume adjustment
- Separate tip ejector
- Clear digital window on the hand grip


## Accessories



### BIOTIPS

- Made from virgin plastic
- Fully autoclavable
- Universal microtips

Cat No.	Description	Pack / Size
BS/003/WS	White (0.2-10 µl)	1x1000 pcs.
BS/003/Y	Yellow (2-200 µl)	1x1000 pcs.
BS/003/B	Blue (200-1000 µl)	1x1000 pcs.
BS/003/WM	White (0.5-5 ml)	1x100 pcs.
BS/003/WL	White (1-10 ml)	1x100 pcs.




### BIOTIPS BOX

- Made from heavy plastic
- The colour coding of the tip may allows immediate tip volume recognition.

Cat No.	Description	Pack / Size
BS/017/W	White (0.2-10 µl)	1x96 pcs.
BS/017/Y	Yellow (2-200 µl)	1x96 pcs.
BS/017/B	Blue (200-1000 µl)	1x96 pcs.

### BIOSTAND



Cat No.	Description	Pack / Size
BS/004/A	BIOSTAND 4holes	1 pcs.
BS/004/B	BIOSTAND 5holes	1 pcs.
BS/004/C	BIOSTAND 2holes	1 pcs.
BS/004/D	BIOSTAND 1holes	1 pcs.

## Order Specification

### FIXED VOLUME

Cat. No.	Volume	Inaccuracy	Imprecision	Tips to be used
BS/013/5	5 µl	± 2.5	< 1.5	Yellow
BS/013/10	10 µl	± 1.2	< 0.8	Yellow
BS/013/20	20 µl	± 1.0	< 0.5	Yellow
BS/013/25	25 µl	± 1.0	< 0.5	Yellow
BS/013/50	50 µl	± 0.8	< 0.4	Yellow
BS/013/100	100 µl	± 0.8	< 0.3	Yellow
BS/013/200	200 µl	± 0.8	< 0.3	Yellow
BS/013/250	250 µl	± 0.8	< 0.3	Blue
BS/013/500	500 µl	± 0.8	< 0.3	Blue
BS/013/1000	1000 µl	± 0.8	< 0.3	Blue
BS/013/5000	5000 µl	± 0.8	< 0.3	White
BS/013/10000	10000 µl	± 0.8	< 0.3	White

### VARIABLE VOLUME

Cat. No.	Volume	Increment	Inaccuracy	Imprecision	Tips to be used
BS/014/200	0.5 - 10	0.1	± 5.0 - 2.0	< 5.0 - 1.0	Yellow
BS/014/210	10 - 100	1.0	± 2.0 - 1.5	< 1.0 - 0.6	Yellow
BS/014/220	5 - 50	1.0	± 1.5 - 2.0	< 0.6 - 0.4	Yellow
BS/014/230	2 - 20	0.1	± 4.0 - 2.0	< 2.0 - 0.8	Yellow
BS/014/240	20 - 200	1.0	± 2.0 - 1.5	< 0.8 - 0.4	Yellow
BS/014/250	100 - 1000	10.0	± 1.5 - 1.2	< 0.6 - 0.4	Blue
BS/014/260	0.5 - 5 ml	100.0	± 5.0 - 2.0	< 0.5 - 1.0	White
BS/014/270	1 - 10 ml	100.0	± 4.0 - 2.0	< 2.0 - 0.8	White

### MULTI-CHANNEL VARIABLE VOLUME

Cat. No.	Volume	Inaccuracy	Imprecision	Tips to be used
BS/014/MV/50	5 - 50	± 2.50 - 0.60	< 1.50 - 0.60	Yellow
BS/014/MV/300	30 - 300	± 4.70 - 1.20	< 1.80 - 0.70	Yellow

NOTE : Special volume pipettes other than specified above may be supplied on demand.

## BIOSYSTEM INTERNATIONAL

Head Office: "India House", 56/6 Popatwadi, 1st Floor, Near Hanuman Temple, Kalbadevi Road, Mumbai 400 002, INDIA. • Telefax: +91-0250-2393979

Work Office: 120, Shakar Industrial Estate, Navghar Road, Vasai (E), Thane 401201, Maharashtra, India. • Telefax: +91-0250-2393979

Email: [info@micropipette.in](mailto:info@micropipette.in) / [biosystem2002@yahoo.com](mailto:biosystem2002@yahoo.com)

Website: [www.micropipette.in](http://www.micropipette.in)