

# Classic Plus

Compare Price Performance Value  
There's none better !



## Every Drop of Precision

- Comfortable finger rest
- Push button with soft rubber
- The weight just 75 gm (Appro.)
- Instant volume setting
- Simple in lab calibration
- Economical price



**BIOSYSTEM**  
'ACCURACY & PRECISION'

CE ISO 9001 : 2000 CERTIFICATE CO.

# Classic Plus

Biosystem newly introduced a (Classic Pro) Model at very affordable price for pathology, medical college, professional educational, laboratory use with the full features of the advance features type micropipettes.



## Additional features of Classic Plus

- Comfortable finger rest
- Push button with soft rubber
- The Weight just 75 gm (Approximately)
- Instant volume setting
- Simple in lab calibration
- Economical price

## Accessories



### BIOTIPS

- Made from virgin plastic
- Fully autoclavable
- Universal microtips

Cat No.	Description	Pack / Size
BS/003/WS	White (0.2-10 µl)	1x1000 pcs.
BS/003/Y	Yellow (2-200 µl)	1x1000 pcs.
BS/003/B	Blue (200-1000 µl)	1x1000 pcs.
BS/003/WM	White (0.5-5 ml)	1x100 pcs.
BS/003/WL	White (1-10 ml)	1x100 pcs.




### BIOTIPS BOX

- Made from heavy acrylic
- Universal pipette stand fitting upto 4 pipette & with tips sprae.

Cat No.	Description	Pack / Size
BS/017/W	White (0.2-10 µl)	1x96 pcs.
BS/017/Y	Yellow (2-200 µl)	1x96 pcs.
BS/017/B	Blue (200-1000 µl)	1x96 pcs.

### BIOSTAND



Cat No.	Description	Pack / Size
BS/004/A	BIOSTAND 4holes	1 pcs.
BS/004/B	BIOSTAND 5holes	1 pcs.
BS/004/C	BIOSTAND 2holes	1 pcs.
BS/004/D	BIOSTAND 1holes	1 pcs.

## Order Specification

### FIXED VOLUME

Cat. No.	Volume	Inaccuracy	Imprecision	Tips to be used
BS/005/5	5 µl	± 2.5	< 1.5	Yellow
BS/005/10	10 µl	± 1.2	< 0.8	Yellow
BS/005/20	20 µl	± 1.0	< 0.5	Yellow
BS/005/25	25 µl	± 1.0	< 0.4	Yellow
BS/005/50	50 µl	± 1.0	< 0.4	Yellow
BS/005/100	100 µl	± 0.8	< 0.3	Yellow
BS/005/200	200 µl	± 0.8	< 0.3	Yellow
BS/005/250	250 µl	± 0.8	< 0.3	Blue
BS/005/500	500 µl	± 0.8	< 0.3	Blue
BS/005/1000	1000 µl	± 0.8	< 0.3	Blue
BS/005/5000	5000 µl	± 0.8	< 0.3	White
BS/005/10000	10000 µl	± 0.8	< 0.3	White

### VARIABLE VOLUME

Cat. No.	Volume	Increment	Inaccuracy	Imprecision	Tips to be used
BS/006/200	0.5 - 10	0.1	± 5.0 - 2.0	< 5.0 - 1.0	White
BS/006/210	10 - 100	1.0	± 2.0 - 1.5	< 1.0 - 0.6	Yellow
BS/006/220	5 - 50	1.0	± 2.5 - 1.4	< 1.5 - 0.6	Yellow
BS/006/230	2 - 20	0.1	± 4.0 - 2.0	< 2.0 - 0.8	White
BS/006/240	20 - 200	1.0	± 2.0 - 1.5	< 0.8 - 0.4	Yellow
BS/006/250	100 - 1000	10.0	± 1.5 - 1.2	< 0.6 - 0.4	Blue
BS/006/260	0.5 - 5 ml	100.0	± 1.5 - 0.5	< 0.6 - 0.4	White
BS/006/270	1 - 10 ml	100.0	± 1.5 - 0.5	< 0.6 - 0.4	White

### MULTI-CHANNEL VARIABLE VOLUME

Cat. No.	Volume	Inaccuracy	Imprecision	Tips to be used
BS/006/MV/50	5 - 50 µl	± 5.00 - 1.50	< 2.00 - 0.80	Yellow
BS/006/MV/300	30 - 300 µl	± 4.70 - 1.20	< 1.80 - 0.70	Yellow

NOTE : Special volume pipettes other than specified above may be supplied on demand.

## BIOSYSTEM INTERNATIONAL

Head Office: "India House", 56/6 Papatwadi, 1st Floor, Near Hanuman Temple, Kalbadevi Road, Mumbai 400 002, INDIA. • Telefax: 0250 2393979  
Email: [info@micropipett.in](mailto:info@micropipett.in) / [biosystem2002@yahoo.com](mailto:biosystem2002@yahoo.com)  
Website: [www.micropipett.in](http://www.micropipett.in)

Work Office: 120, Shakar Industrial Estate, Navghar Road, Vasai (E), Thane 401201, Maharashtra, India. • Telefax: 0250 2393979